





REV	ECO	DESCRIPTION	DATE	DFTM	CHKD
E	13949	CORRECTED OVERALL HEIGHT DIMENSION	1/23/2019	TH	BK



MODEL DIMS			
MODEL	DIM "A"	DIM "B"	HEPA FILTER SIZE
NU-340-630	32 [813]	76 3/8 [1940]	70 [1778] X 30 [762] X 3 [76]

* OVERALL DIMENSIONAL TOLERANCE ±1/4" [6.35]
ALL OTHER DIMENSIONS ±1/8" [3.17]

TITLE NU-340-600 HORIZONTAL AIRFLOW CONSOLE					
DFTM LS	DATE 2/8/2018	CHKD BP	SHEET 1	OF 1	
DRAWING NUMBER BCD-19667			E		

SECTION A-A

INCH
MILLIMETER