



## ACID STORAGE CABINETS

NuAire Acid Storage Cabinets are constructed from all stress-relieved, fully seam-welded, white corrosive-resistant polypropylene.

The outer cabinet shelf, kick plate and adjustable shelves with drip-proof lips are constructed from ½ inch polypropylene.

Each cabinet incorporates easy-to-clean shelves, handles, hinges, and screws also made of 100% polypropylene.

Adjustable slide front air vents and removable top air vent plugs assure acid storage cabinets have adequate ventilation for safety.

DOOR LOCKS AND CASTERS ARE OPTIONAL.

## ACID STORAGE CABINETS



MODEL NUMBER	SIZE (W x H x D)	PRICE
NU-53-3660	36" x 60" x 18"	\$1,915.00
NU-53-4260	42" x 60" x 18"	1,975.00
NU-53-4860	48" x 60" x 18"	2,035.00
NU-53-3684	36" x 84" x 22"	\$2,725.00
NU-53-4284	42" x 84" x 22"	2,965.00
NU-53-4884	48" x 84" x 22"	3,200.00

**OPTIONS:** Casters, Vent Location, Key Locks

## POLYPROPYLENE TABLES



MODEL NUMBER	SIZE (W x H x D)	PRICE
NU-98-3636	36" x 36" x 24"	\$1,250.00
NU-98-4836	48" x 36" x 24"	1,470.00
NU-98-6036	60" x 36" x 24"	1,600.00
NU-98-7236	72" x 36" x 24"	1,795.00

**OPTIONS:** Casters, Epoxy Resin Tabletop

## PEG BOARDS



### WHITE/BLACK W/PEGS & 3" DRAIN TROUGH




MODEL NUMBER	SIZE (W x H)	PRICE
NU-56-3030	30" x 30"	\$450.00
NU-56-3036	30" x 36"	\$500.00

## WITH NUAIRE YOU KNOW



[www.nuaire.com](http://www.nuaire.com)

2100 Fernbrook Lane  
Plymouth, MN 55447  
U.S.A.  
Phone: 763.553.1270  
Fax: 763.553.0459  
1.800.328.3352



PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE. U.S. DOLLARS ONLY. SUBJECT TO APPLICABLE TARIFFS AND TAXES. PRICES EFFECTIVE JANUARY 2006.