

**TECHNICAL BULLETIN  
ACCESSORY INFORMATION**




2100 Fernbrook Lane  
Plymouth, MN 55447-4722  
U.S.A.

Phone: 763.553.1270  
Fax: 763.553.0459  
Toll Free: 1.800.328.3352

**NU-976 Disposable Forearm Pads  
Installation Instructions**

The enclosed disposable forearm pads are made from 1/4" (6mm) thick closed cell cross-link polyethylene material. The disposable pads were specifically designed and tested for NuAire Biological Safety Cabinets. Please consult original equipment manufacturer for use on other than NuAire Biological Safety Cabinets. The pad attachment is done by peeling off the strip tape release paper, and placing pad as illustrated below. The disposable pads are intended for single use only and should be treated as biohazard waste.



Additional pads can be ordered by calling NuAire Small Part Order Entry @ 1-800-328-3352.