

**TECHNICAL BULLETIN
ACCESSORY INFORMATION**



2100 Fernbrook Lane
Plymouth, MN 55447-4722
U.S.A.

Phone: 763.553.1270
Fax: 763.553.0459
Toll Free: 1.800.328.3352

**NU-977-001 Rotating Table
Installation Instructions**

The enclosed rotating table or Lazy-Susan table was specifically designed for use within a Biological Safety Cabinet (BSC). The table is 12 inches (305mm) in diameter and 1/2 inch (13mm) high. The table is constructed from stainless steel with a galvanized bearing system capable of handling 100 lbs. (45 Kg).

The rotating table is autoclavable or can be disinfected with a decontaminate such as alcohol
**DO NOT USE CHLORINATED OR HALOGEN MATERIALS, AS THEY MAY
DAMAGE THE STAINLESS STEEL.**