





REV	ECO	DESCRIPTION	DATE	DFTM	CHKD
D	13949	CORRECTED OVERALL HEIGHT DIMENSIONS	1/23/2019	TH	BK



MODEL DIMS			
MODEL	DIM "A"	DIM "B"	HEPA FILTER SIZE
NU-340-430	32 [813]	76 3/8 [1940]	46 [1168] X 30 [762] X 3 [76]
NU-340-436	38 [965]	82 3/8 [2092]	46 [1168] X 36 [914] X 3 [76]

* OVERALL DIMENSIONAL TOLERANCE ±1/4" [6.35]
ALL OTHER DIMENSIONS ±1/8" [3.17]

TITLE: NU-340-400 HORIZONTAL AIRFLOW CONSOLE

DFTM	DATE	CHKD	SHEET	OF
LS	3/14/2018	BP	1	1

DRAWING NUMBER: BCD-19803

D

INCH
MILLIMETER