





REV	ECO	DESCRIPTION	DATE	DFTM	CHKD
A	13949	RELEASED TO PRODUCTION	1/23/2019	TH	BK



MODEL DIMS			
MODEL	DIM "A"	DIM "B"	HEPA FILTER SIZE
NU-340-530	32 [813]	80 1/2 [2045]	58 [1473] X 30 [762] X 3 [76]
NU-340-536	38 [965]	86 1/2 [2197]	58 [1473] X 36 [914] X 3 [76]

* OVERALL DIMENSIONAL TOLERANCE ±1/4" [6.35]
ALL OTHER DIMENSIONS ±1/8" [3.17]

TITLE: NU-340-500 HORIZONTAL AIRFLOW CONSOLE W/ CASTORS

DFTM TH	DATE 1/23/2019	CHKD BK	SHEET 1 OF 1
DRAWING NUMBER		BCD-20629	

A

SECTION A-A

INCH MILLIMETER