




Parts List		
ITEM	QTY	DESCRIPTION
1	2	WELDMENT, NU-813/819 (S-20) BASE STAND END 36"
2	1	NU-813/819-600 (S-20) BASESTAND CROSS BRACE
3	2	Lateral Brace
4	2	Nutsert 1/4-20
5	4	1/2-13 2-1/2" Leg Leveler
6	2	Bolt G8 3/8-16 X 3 Nickel
7	2	Washer 3/8 SAE Nickel
8	2	Lock Washer 3/8 Internal Tooth SST
9	4	Bolt G2 1/4-20 X 2-1/2 SST
10	4	Acorn Nut 1/4-20 (Low Crown) SST
11	6	Lock Washer 1/4 Internal Tooth SST
12	10	Washer 1/4 SST
13	2	Bolt G2 1/4-20 X 3/4 SST

REV	ECO	DESCRIPTION	DATE	DFTM	CHKD
A	000925	RELEASED	1/31/2022	RN	DF

1/4-20 NUTSERT CENTER HOLE
BOTH BASE STAND END WELDMENTS



		PROPRIETARY THE INFORMATION CONTAINED HEREIN IS THE EXCLUSIVE PROPERTY OF NUAIRE INC. AND IS NOT TO BE DIVULGED OR USED IN ANY MANNER WITHOUT THE EXPRESS WRITTEN PERMISSION OF NUAIRE INC.	
		DATE	1/31/2022
DFTM	RN		
CHKD	DF		
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES -TOLERANCES- DECIMALS ± .032 ANGLES ± 2°		MATERIAL	AS NOTED
		PART NUMBER	CD-000136
		PROJECT NUMBER	
			SHEET 1 OF 1