



POLYPROPYLENE

RUST FREE ACID STORAGE CABINETS



NuAire Acid Storage Cabinets are constructed from all stress-relieved fully seam-welded, white corrosive-resistant polypropylene.

The outer cabinet shelf, kick plate, and adjustable shelves with drip-proof lips are constructed from $\frac{1}{4}$ inch polypropylene.


Each cabinet incorporates easy-to-clean shelves, handles, hinges, and screws are also made of 100% polypropylene.

Adjustable slide front air vents and removable top air vent plugs assure acid storage cabinets have adequate ventilation for safety.



Best Products. Best Performance. Best Protection.

Polypropylene Acid Storage Cabinets



Specifications

Overall Dimensions	Width (A)	Height (B)	Depth (C)
NU-53-3660	36" (914 mm)	60" (1524 mm)	18" (457 mm)
NU-53-4260	42" (1067 mm)	60" (1524 mm)	18" (457 mm)
NU-53-4860	48" (1219 mm)	60" (1524 mm)	18" (457 mm)
NU-53-3684	36" (914 mm)	84" (2134 mm)	22" (559 mm)
NU-53-4284	42" (1067 mm)	84" (2134 mm)	22" (559 mm)
NU-53-4884	48" (1219 mm)	84" (2134 mm)	22" (559 mm)

Options: Casters, Vent Location, Key Locks, Custom Configurations

Polypropylene Laboratory Tables & Carts



NU-98: Polypropylene Table



NU-96: Laboratory Cart

Specifications

Overall Dimensions	Width	Height	Depth
NU-98-3636	36" (914 mm)	36" (914 mm)	24" (610 mm)
NU-98-4836	48" (1219 mm)	36" (914 mm)	24" (610 mm)
NU-98-6036	60" (1524 mm)	36" (914 mm)	24" (610 mm)
NU-98-7236	72" (1829 mm)	36" (914 mm)	24" (610 mm)
NU-96	25" (635 mm)	26" (660 mm)	19 3/4" (502 mm)

Polypropylene Peg Boards



Specifications

Overall Dimensions	Width	Height
NU-56-3030	30" (762 mm)	30" (762 mm)
NU-56-3036	30" (762 mm)	36" (914 mm)

